

Assays Using MCF-7 Cells

Reference	Deckers et al. (2000)	Schoonen et al. (1995)
Characteristics of Cells		
<i>Cell line</i>	MCF-7	MCF-7
<i>Source of cell line</i>	human mammary tumor	human mammary tumor
<i>Whole cells/cytosol</i>	cytosol	cytosol
Preparation of Cells for Assay		
<i>Serum source</i>	fetal calf serum	fetal calf serum
<i>Serum stripping method</i>	charcoal treated serum	charcoal treated serum
<i>Residual androgen in serum</i>	n.p.	n.p.
<i>No. treated cells/No. or weight of cells homogenized</i>	1 gm cells	1 gm cells
<i>Treatment vessel used</i>	n.p.	n.p.
<i>Preparation of cell homogenate</i>		
<i>volume</i>	5 ml	5 ml
<i>buffer</i>	TrisHCl pH 7.4 + EDTA, dithioerythritol, molybdate	TrisHCl pH 7.4 + EDTA, dithioerythritol, molybdate
<i>method</i>	Dounce homogenizer	Dounce homogenizer
<i>time; temperature</i>	n.p.	n.p.
<i>Centrifugation of homogenate</i>	1,000,000N/kg	1,000,000N/kg
<i>Protein concentration of cytosol</i>	n.p.	n.p.
<i>Storage</i>	n.p.	n.p.
<i>Final protein concentration</i>	n.p.	n.p.
Competitive binding assay		
<i>Reference ligand</i>	5 α-Dihydrotestosterone	5 α-Dihydrotestosterone
<i>Volume and concentration of reference ligand</i>	1.9 nM	1.9 nM
<i>Specific activity of labelled reference ligand</i>	5.3 TBq/mmol	4070 GBq/mmol
<i>Volume and concentration of cold ligand</i>	n.p.	n.p.
<i>Final concentration of reference ligand</i>	1.9 nM	1.9 nM
<i>Volume of competing ligand</i>	n.p.	n.p.
<i>Concentration range of competing ligand</i>	0.1 - 10000 nM	0.1 - 10000 nM
<i>Volume of cytosol</i>	1:5 dilution	1:5 dilution
<i>Volume of buffer</i>	n.p.	n.p.
<i>Type of buffer used</i>	n.p.	n.p.
<i>Replicates</i>	6 or more	2
<i>Time of incubation</i>	overnight	overnight
<i>Temperature of incubation</i>	4° C	4° C

Assays Using MCF-7 Cells

Reference	Deckers et al. (2000)	Schoonen et al. (1995)
Separation of ligand		
Volume and type of slurry	dextran-charcoal	dextran-charcoal
Buffer for slurry	TrisHCl pH 7.4 + EDTA, dithioerythritol, molybdate	TrisHCl pH 7.4 + EDTA, dithioerythritol, molybdate
Incubation time and temperature	10 min, 4° C	10 min, 4° C
Time of vortexing	n.a.	n.a.
Centrifugation speed	8000N/kg	8000N/kg
Centrifugation time and temperature	5 min	5 min
Resuspension volume and buffer for pellet	n.p.	n.p.
No. of washes	n.p.	n.p.
Extraction of label	centrifugation	centrifugation
Incubation time and temperature	n.p.	n.p.
Vortexing during incubation time	n.p.	n.p.
Centrifugation time and temperature	n.p.	n.p.
Volume added for reading	n.p.	n.p.
Volume of fluor	n.p.	n.p.
Type of fluor	n.p.	n.p.
Instrumentation	Topcount microplate scintillation counter	n.p.
Measurement	n.p.	n.p.
Blank without competitor	n.p.	n.p.
Reading of blank	n.p.	n.p.
Blank subtracted?	n.p.	n.p.
Range of standard curve of reference ligand	n.p.	n.p.
Nonspecific binding measured?	yes	n.p.
Subtraction of nonspecific binding	yes	n.p.
Data calculations		
Data plotted as	n.p.	n.p.
Data calculated	specific binding	IC ₅₀
Calculation of RBA	yes	yes
Test substances		
Solvent used	ethanol	ethanol
No. of samples/ dose	n.p.	n.p.
No. of times assay repeated	from 6 to 34	2
Abbreviations: n.a. = not applicable; n.p. = not provided; RBA = relative binding affinity		